



MACK'S®
audibles®
Personal Sound Amplifier

**Amplify Sound.
Enhance Life!**

What is Mack's® Audibles®?

- Mack's® Audibles® is a wearable electronic product (PSAP*) designed to amplify sound for those who may need an occasional volume "boost" to hear more clearly.

**a PSAP is a Personal Sound Amplifier - it is not a hearing aid.*

- Mack's® Audibles® are unique in that the sound amplifier is positioned in the ear canal, providing optimal volume and sound clarity where it matters most – in the ear.
- Clean, simple, fashionable design.
- Comfortable to wear. Easy to use.
- If your patients sometimes need to amplify the world around them in business, entertainment, educational or social situations, Mack's® Audibles® are a simple solution to help them hear what they've been missing.



Mack's® Audibles® Personal Sound Amplifier Technical Specs



Digital and analog technology supplies clear, crisp sound enhancement

**Maximum output = 116 dB
(limits loudest sounds)**

Powerful sound gain = up to 43 dB

Total harmonic distortion = <3 dB

Two listening modes for louder or quieter environments

Fits either ear

MACK'S®
audibles®
Personal Sound Amplifier